



Teeth Whitening Premium X3 Kit

Art no: 121897

Up to 6 shades whiter teeth in 14 days

Claims:

- Up to 6 shades whiter teeth clinically proven
- Safe and without side effects
- Patented Gel formula with Dual Boost^{x2} & ^{x3} Technology that cleans with detergents, whitens with sodium bicarbonate and coco-sulphate at the same time + X3 gel with 0,1% peroxide for an even deeper whitening effect.
- Kit includes Tooth Gloss that protects against new stains and gives better shine to your teeth.
- Customized and reusable mouth trays
- Up to 14 days treatment

Product description:

The X3 Premium Kit is a complete system to whiten your teeth at home fast, easy and gentle. It contains 8 ml Peroxide Free Dual Boost Gel^{x2} + 8 ml 0,1% Peroxide Dual Boost^{x3} gel, which is used together with the thermoform mouth trays that are individually fitted using hot water.

The product concept is to clean and whiten the teeth gently with the Dual Boost^{x2} gel (Peroxide Free) for 3-4 days and then use the Dual Boost^{x3} gel (0,1% Peroxide) for another week or until desired result is achieved.

The Premium Kit also contains a Tooth Gloss which is used after the teeth whitening for protection of the enamel and giving the teeth an extra gloss.

Working mechanism in short:

The Dual Boost^{x2} gel.

The detergent technology creates negatively charged ionic foaming bubbles that work via ionic attraction to lift stain molecules away from teeth without increasing tooth sensitivity associated with typical bleaching gels that dentists use. The product also has special sodium bicarbonate whiteners, solvents and wetting agents for a deeper second wave of whitening activity that begins to work after 8 minutes. The detergent and whiteners together gives a Dual Boost activity for a more effective result.

The Dual Boost^{x3} gel.

The X3 gel is basically the same as X2 but with added 0,1% hydrogen peroxide. The peroxide concentration is within EU legal limit and low enough not to cause any side effects. Adding peroxide will improve the whitening effect with approximately one shade point. The oxidizing effect will further improve the cleaning effect and will make the teeth appear whiter.

Unique Competitive Advantages Dual Boost^{x2} and Dual Boost^{x3} gel:

Group Beconfident Europe AB

Tennisvägen 31, SE-26933 Båstad, Sweden

E-mail: info@beconfident.info Phone: +46-(0)702 56 50 00 VAT no: SE 556853-0041-01

Corporate site: www.beconfident.info, www.beconfident.co.uk, www.beconfident.se & many more

Beconfident- Top Quality Teeth whitening and Skincare systems are designed in Sweden & Made in USA and Sweden. Beconfident® + WhiteAmin® + SimpleSmile® are registered trademarks in Europe and Beconfident® is a registered trademark in United States of America. IOS and Android Intelligent Apps belong to Beconfident Skincare AB. Other brands are pending.





- Unique Dual Boost^{x2} technology for a safe process and efficient result that uses 2 different mechanisms to whiten teeth – Cleaning and Whitening.
- Dual Boost^{x3} with added 0,1% peroxide will further improve the cleaning and whitening effect. Fast Dual Whitening Activity -> 1-2 shades in 30 minutes up to 5-6 (After full treatment)
- Whitening treatments last longer because stains are actually removed.
- Both gel formulations made of eco healthy safe ingredients derived from Coconuts.
- Three mouth trays that are individually fitted to optimize the contact between the gel and tooth surface. One tray extra in reserve.
- Non-Abrasive. Does not scratch teeth enamel and works safely for all users.
- Safe legal concentration of peroxide for improved effect.
- Includes free Teeth Whitening app to guide and support customer during treatment in order to achieve the best result.

Unique Competitive Advantages Tooth gloss gel:

- Tooth gloss protects enamel from discoloration after the whitening process
- Brightens teeth (higher luster shine)
- Delicious flavor that refreshes breath
- Easy to carry along in a hygienic container

Ingredients:

Syringe 1 Dual Boost^{x2} Gel:

Glycerin, Sodium Bicarbonate, Sodium Coco-Sulfate, Mentha Piperita Extract, Xylitol, Chondrus Crispus, Sorbitol, Aloe Leaf Juice, Pomegranate Seed Extract, Chamomile Flower Extract.

Syringe 2 Dual Boost^{x3} Gel:

Glycerin, Sodium Bicarbonate, Sodium Coco-Sulfate, Mentha Piperita Extract, Xylitol, Chondrus Crispus, Sorbitol, Aloe Leaf Juice, Pomegranate Seed Extract, Chamomile Flower Extract, Carbamide Peroxide.

Tooth Gloss Mint:

Aqua, Poloxamer 407, Alcohol, Carbopol, Sucralose, Sodium Saccharin, L-Menthol, Sodium Benzoate, Aroma.

The Premium kit contains the following items:

- 1) 1 x 8 ml Syringes with Dual Boost^{x2} Peroxide Free Gel
- 2) 1 x 8ml Syringe with Dual Boost^{x3} 0,1% Peroxide Gel
- 3) 1 x 10 ml Tooth Gloss
- 4) 2 shapeable mouth trays
- 5) Tooth Gloss applicator pen
- 6) Instruction sheet and a shade guide.

Beconfident is not conducting any testing on animals.

Product Data:

Art No	Dimension HxWxD	Gross Weight	Net Weight	Units Display Box	Units mastercart.	Unit EAN
--------	--------------------	--------------	------------	----------------------	----------------------	----------

Group Beconfident Europe AB

Tennisvägen 31, SE-26933 Båstad, Sweden

E-mail: info@beconfident.info Phone: +46-(0)702 56 50 00 VAT no: SE 556853-0041-01

Corporate site: www.beconfident.info, www.beconfident.co.uk, www.beconfident.se & many more

Beconfident- Top Quality Teeth whitening and Skincare systems are designed in Sweden & Made in USA and Sweden. Beconfident® + WhiteAmin® + SimpleSmile® are registered trademarks in Europe and Beconfident® is a registered trademark in United States of America. IOS and Android Intelligent Apps belong to Beconfident Skincare AB. Other brands are pending.

beconfident[®]

Product Data Sheet



	(mm)					
1218xx	175x120x33	153g	121g	4	102	See below
CPNP Code	GMP ISO 22716	Country of Origin	Quality Test Pharmacies	Box bio- degradable	Pallet	RRP Euro incl VAT(+/-)
	Follows	USA		Yes	12	79,90 Euro

Article numbers

EU NORTH -51 (EN, SV, NO/DK, FI)
EU CENTRAL -52 (EN, FR, NL, DE, IT)
EU EAST -53 (EN, CZ, SK, PL)

Group Beconfident Europe AB

Tennisvägen 31, SE-26933 Båstad, Sweden

E-mail: info@beconfident.info Phone: +46-(0)702 56 50 00 VAT no: SE 556853-0041-01

Corporate site: www.beconfident.info, www.beconfident.co.uk, www.beconfident.se & many more

Beconfident- Top Quality Teeth whitening and Skincare systems are designed in Sweden & Made in USA and Sweden. Beconfident[®] + WhiteAmin[®] + SimpleSmile[®] are registered trademarks in Europe and Beconfident[®] is a registered trademark in United States of America. IOS and Android Intelligent Apps belong to Beconfident Skincare AB. Other brands are pending.